

# ARIANE

## PUSH-IT PACK : 2 FORMULATIONS TO MIX



- PP / PE / ELASTOMER

Ø20 - Ariane Push-it



15 ml  
H: 147 mm



### ZOOM ON THE PRODUCT

This system allows' to the final consumer to activate the formulation (mixing both products)

Possibility to have a nomad pack with a screwing closure

Packaging delivered in 3 parts (a ampoule with integrated applicator tip, a tank and a valve)

For aqueous formulations from 1 to 3000 centipoises




Available with different applicator tips adapted to your formulas

A - Integrated applicator tip  
B - Headwear applicator tip

The Ariane Pack allows two formulations to be activated at the time of use:

- 1 REMOVE tamper tab
- 2 PRESS the cap of the tank for release the first formula: powder or liquid
- 3 MIX the two formulas
- 4 BREAK the head of the ampoule
- 5 APPLY pressure on the cap to distribute the product: drop by drop or larger dose



### DECORATION

Screen (silk) printing  
Hot stamping

You can create your own artwork without technical limit.

Ø20 ARIANE PUSH-IT  
DISTRIBUTION : Liquids, Powders, Serum, Cream